










The Design					
The Product Line	The Blue Line				
One-Knob Data Entry or Comfortable Touchpad					
					
Diameter of the Chamber in mm	254	254	355	355	355
Model	MicroCen M	MicroCen MR	UniCen M	UniCen MR	UniCen HR
Table / Floor	Table top	Table top	Table top	Table top	Table top
Max. Speed	16.000	22.000	18.000	24.000	30.000
g-Force	20.464	37.878	25.356	45.078	70.484
Refrigerated	no	yes	no	yes	yes
Max. Volume Swing-out in ml	24 x 2,2 ml	24 x 2,2 ml	4 x 200	4 x 200	4 x 200
Max. Volume Fixed-Angle in ml	12 x 10	12 x 10	6 x 100	6 x 100	6 x 100
Max. Capacity 1,5/2,2 ml Tubes	30	30	48	48	48
Max. Capacity 15 ml Falcon® Tubes	n.a.	n.a.	30	30	30
Max. Capacity 50 ml Falcon® Tubes	n.a.	n.a.	8	8	8

Herolab reserves the right to change technical specification without prior notice. 04/2016

